



SHANGHAI JOY CROWN INDUSTRY CO., LTD.

Always Your Best Choice

[HTTP://WWW.JOYCROWN.COM](http://www.joycrown.com)





Brief of Company:

Shanghai Joy Crown Industry Co., Ltd is specializing in production of medical gauze products, elastic bandage series, non woven care products, surgical plaster, wound dressing, first aid kits and other medical disposables. Our company is in Shanghai, the busiest port in the China Mainland. We have almost 20 years' experience in medical supplies. We provided many kinds of medical products, Class I and Class II to Europe, America, South-America, Middle-East, East-South Asia and Africa in OEM. And we we are capable of developing new products for customers as they required.

Our products are up to the standard of BP 2001 and USP Standard. The main products are including all kinds of medical gauze products, absorbent zigzag gauze, non-sterile and sterile gauze rolls, gauze swabs, lap sponges, gauze bandages, gauze balls, and first aid bandages, elastic bandages series, pure (100%) cotton stretched bandages, and crepe bandages made by 67% cotton and 34% wool, high elastic bandages, elastic adhesive bandages, zinc oxide adhesive plasters and surgical tapes (non woven tape, PE tape, silk tape). PU Incise Drapes and wound dressings, P.O.P bandage and slabs, casting tapes, casting paddings, and all kinds of non-woven care products (non woven sponges, face masks, caps, shoe covers, isolation gowns, surgical gowns and so on).

Also, we are dealing with foam cervical collars, adjustable cervical collars, hard plastic and polyethylene outer shell, dry wipes, wet wipes, surgical scissors, thermometers & digital thermometers, foam sponges, first aid kits and surgical kits, alcohol swabs, mouth to mouth masks (rescuer masks), emergency blankets, toiletry bags, pregnancy test kits, syringes, infusion sets.

Shanghai Joy Crown Industry Co., Ltd. is moving aggressively by continuing its commitment to quality products, competitive prices, fast delivery time and effective customer service. Joy Crown is always your best choice in business!

Our Reference:

- 3M
- Smith & Nephew
- Clinisupplies UK
- Roval Med S.R.L



Always Your Best Choice



SHANGHAI JOY CROWN INDUSTRY CO., LTD.

CE 0197



Always Your Best Choice



GAUZE SERIES

<http://www.joycrown.com>



Absorbent Gauze Zigzag
 Ref. No.: JC01001
 40's, 32's, 21's; 13, 17, 20 threads
 43x27, 30x30, 28x24, 24x20,
 26x17, 19x15, 15x18, 20x12
 Size: 36" x 40m, 36" x 100m(yds)
 48" x 100m(yds); 80cm x 100m;
 100cm x 100m
 Packing: zigzag folded, poly bag package



Gauze Roll
 Ref. No.: JC01002
 40's, 32's, 21's; 13, 17, 20 threads
 43x27, 30x30, 28x24, 24x20,
 26x17, 19x15, 15x18, 20x12
 Size: Width: 36", 48", 100cm
 Jumbo roll: 2000/3000m flat roll
 Medium roll: 100m(yds) flat roll 4ply



Gauze Swab
 Ref. No.: JC01003
 40's, 32's, 21's; 13, 17, 20 threads
 43x27, 30x30, 28x24, 24x20,
 26x17, 19x15, 15x18, 20x12
 Size: 5x5cm, 7.5x7.5cm,
 10x10cm, 10x20cm
 Ply: 8ply, 12ply, 16ply
 * with or without X-ray
 * non-sterile



Sterile Gauze Swab
 Ref. No.: JC0100391
 40's, 32's, 21's; 13, 17, 20 threads
 43x27, 30x30, 28x24, 24x20,
 26x17, 19x15, 15x18, 20x12
 Size: 5x5cm, 7.5x7.5cm,
 10x10cm, 10x20cm
 Ply: 8ply, 12ply, 16ply
 Package: 1pc, 2pcs, 5pcs, 10pcs/pack
 * with or without X-ray
 * sterile



Lap Sponge(Abdominal Pad)
 Ref. No.: JC01004
 40's, 32's, 21's; 13, 17, 20 threads
 43x27, 30x30, 28x24, 24x20,
 26x17, 19x15, 15x18, 20x12
 Size: 30cm x 30cm; 45cm x 45cm;
 4ply, 6ply, 8ply, 12ply.
 *with or without X-ray
 *cotton loop: white, blue, green
 *Multiples, sewn
 *non-sterile



Sterile Lap Sponge(Abdominal Pad)
 Ref. No.: JC0100401
 40's, 32's, 21's; 13, 17, 20 threads
 43x27, 30x30, 28x24, 24x20,
 26x17, 19x15, 15x18, 20x12
 Size: 30cm x 30cm; 45cm x 45cm;
 4ply, 6ply, 8ply, 12ply.
 *with or without X-ray
 *cotton loop: white, blue, green
 *multiples, sewn
 *sterile; 5pcs/pouch; double package



Pre-washed Lap Sponge
 Ref. No.: JC0100402
 40's, 32's, 21's; 13, 17, 20 threads
 43x27, 30x30, 28x24, 24x20,
 26x17, 19x15, 15x18, 20x12
 Size: 30cm x 30cm; 45cm x 45cm;
 4ply, 6ply, 8ply, 12ply.
 *sterile or non-sterile
 *with or without X-ray
 *cotton loop: white, blue, green
 *Multiples, sewn



Gauze Bandage
 Ref. No.: JC01005
 Material: 13 threads, 17 threads,
 20 threads, 40's, 32's, 21's
 28x24, 26x18, 19x15
 Length: 4m, 5m, 10m
 Width: 5cm, 7.5cm, 10cm, 12.5cm,
 15cm, 20cm, 30cm, 120cm
 Packing: kraft paper



Sterile Gauze Bandage
 Ref. No.: JC0100501
 Material: 13 threads, 17 threads,
 40's, 32's, 21's
 28x24, 26x18, 19x15
 Length: 2.7m, 4m, 5m, 10m
 Width: 5cm, 7.5cm, 8cm, 10cm,
 12.5cm, 15cm, 20cm,
 30cm, 120cm
 Packing: individual sterile pack



Gauze Ball
 Ref. No.: JC01007
 Material: 40's, 30x20; 24x20,
 20x18, 19x15
 Size: 10x10cm, 15cmx15cm, 20x20cm,
 30x30cm, 40x40cm
 * sterile or non-sterile
 * with or without X-ray
 Packing: non-sterile - 100pcs/PP bag
 sterile - 5pcs/sterile bag



GAUZE & COTTON

<http://www.joycrown.com>



**Krinkle Gauze Roll
(Fluff Roll, Gauze Bandage Roll)**

Ref. No.: JCO1009
Material: 100% cotton, absorbent gauze bandage
Size: 4.5'x144'(4yds, stretched)
*Materials are sterilized and packed in accordance with FDA guidelines
Packing: individual pouch pack



Absorbent Cotton Roll

Ref. No.: JCO2001
Material: 100% cotton with or without interleave paper
Size: 25g, 50g, 100g, 150g, 200g, 250g, 500g, 1000g
Packing: individual pack with kraft paper
*with or without interleave paper



Triangular Bandage

Ref. No.: JCO1010
Material: 100% cotton, unbleached
40's, 60x68, 90x48, 43x27
Size: 90cmx90cmx127cm
95cmx95cmx136cm
100cmx100cmx141cm
106cmx106cmx150cm
Packing: individually packed with 2 safety pins



Cotton Ball

Ref. No.: JCO2002
Material: 100% cotton
Size: 0.5g/pc, 1.0g/pc
Packing: 500pc/bag, 5pc/sterile pack
*sterile or non-sterile



Green Gauze Swab

Ref. No.: JCO1011
40's, 32's, 21's: 13,17,20threads
43x27, 30x30, 28x24, 24x20,
26x17, 19x15, 18x15, 20x12
Size: 5x5cm, 7.5x7.5cm,
10x10cm, 10x20cm
Ply: 8ply, 12ply, 16ply
*sterile or non-sterile



Dental Cotton Roll

Ref. No.: JCO2003
Specially designed for dental use.
Material: 100% cotton
Size: 6x38mm, 10x38mm,
12x38mm
Packing: 50 rolls/bundle



**Burn Sheets
(Cotton Pads)**

Ref. No.: JCO1012
Especially for burning care
20cm extra gauze on both sides
with 1cm thick cotton inside
Size: sheet - 25cmx50cm,
cotton - 25cmx25cmx1cm
Packing: 2pc/sterile pack
*sterile



Cotton Tips

Ref. No.: JCO2004
Material: Cotton tips with stick
Size: 7.5cm, 15cm
*single / double tip
*wood stick / plastic stick
Packing: 100 pcs/bag,
100pc/plastic can



**Paraffin Gauze Dressing
(Vaseline Gauze)**

Ref. No.: JCO1013
Material: gauze coated with paraffin jelly
The paraffin gauze dressing is non-adherent,
non allergenic gamma sterilized, which helps
in the speedy recovery of wounds.
Used in the treatment of ulcers, burns, skin
grafts and various traumatic injuries.
Size: 10cmx10cm, 10cmx20cm, 10cmx30cm,
10cmx40cm, 10cmx7m
Packing: individual pack, sterile



Eye Pad

Ref. No.: JCO2005
Oval-shaped eye pads made of
absorbent cotton filler and enclosed in
a non-woven fabric.
Suitable for eye injury care.
Size: 5.5cm x 7.5cm, 6cm x 8cm
Packing: 1pc/pack, 2pc/pack
*sterile



NON-WOVEN SERIES

<http://www.joycrown.com>



Non-woven Swabs

Ref. No.: JC03001

Material: spun lace non-woven fabric
70% viscose, 30% polyester
Weight: 30gsm, 40gsm, 50gsm
Main size: 5cmx5cm, 7.5cmx7.5cm,
10cmx10cm, 10cmx20cm
4plys, 6plys
Packing: 100pcs/paper pack



Non-woven Triangular Bandage

Ref. No.: JC03004

Material: non-woven material,
Weight: 25gsm, 30gsm, 50gsm
Sizes: 90cm x 90cm x 127cm
100cmx100cmx141cm
Packing: individually packed with
2 safety pins.



Sterile Non-woven Swabs

Ref. No.: JC0300101

Material: spun lace non-woven fabric
70% viscose, 30% polyester
Weight: 30gsm, 40gsm, 50gsm
Main size: 5cmx5cm, 7.5cmx7.5cm,
10cmx10cm, 10cmx20cm
4plys, 6plys
Packing: sterile pack, 1pc/pkt.
2pcs/pkt., 5pcs/pkt.



Dry Wipe

Ref. No.: JC03005

Material: non-woven
Weight: 40gsm
Size: 20cmx20cm-1ply
30cmx18cm-1ply
Packing: 30 pcs/polybag
200bags/carton



Non-woven Swabs with X-ray

Ref. No.: JC0300102

Material: spun lace non-woven fabric
70% viscose, 30% polyester
Weight: 30gsm, 40gsm, 50gsm
Main size: 5cmx5cm, 7.5cmx7.5cm,
10cmx10cm, 10cmx20cm
4plys, 6plys
With X-ray detectable
Packing: 100pcs/paper pack



Combine Pad

Ref. No.: JC03006

Highly absorbent dressing.
Used for absorbing blood and keeping
wounds clean.
Material: non-woven facing with
cellulose or cotton filler.
Size: 5"x9", 8"x10, 10cmx10cm,
10cmx20cm, 20cmx20cm, 20cmx40cm
* non sterile



Non-woven Lap Sponge

Ref. No.: JC03002

Material: spun lace non-woven fabric
70% viscose, 30% polyester
Weight: 30gsm, 40gsm, 50gsm
Size: 10cm x 60cm-4ply
30cm x 30cm-4ply
*with or without X-ray



Under Pad

Ref. No.: JC03021

For incontinent care, medical use, maternity
care, or baby care
Material: Nonwoven + paper tissue + full
pulp + SAP + paper tissue + PE film
Size: Large: 90x60cm,
Medium: 60x60cm,
Small: 60x45cm
Packing: 15pcs/bag



Non-woven Ball

Ref. No.: JC03003

Material: spunlace non-woven fabric
(70% viscose, 30% polyester)
Weight: 30gsm, 40gsm
Sizes: 10cmx10cm, 20cmx20cm,
25cmx25cm, 30cmx30cm
*with or without X-ray
*sterile or non sterile



Non-Adherent Wound Dressing

Ref. No.: J03008

Protect the wound site while promoting
the healing process. The non-stick
facing allows the dressing to absorb
but not adhere to the wound. The easy
removal of the pad minimizes the
trauma to the wound site.
Size: 2"x3", 3"x4", 3"x8"
Packing: Individual sterile pack



NON-WOVEN SERIES

<http://www.joycrown.com>



Alcohol Pad

Ref. No.: JC03009
Components: 70% Isopropyl alcohol
Disinfection, antibacterial, good
for hygienic and sanitary use,
suitable use at home, hospital,
public.
Material: non-woven fabric
Size: 50mmx50mm, 30mmx60mm
Packing: individual sterile pack



Non-woven Face Mask

Ref. No.: JC03013
Material:
2ply: 18g PP+27g PP
3ply: 18g PP+25g filter+27g PP
Size: 17.5cmX9.5cm
*with earloop or tie
Packaging: 50 pcs/box
Color: green, white, blue



Non-woven Isolation Gown

Ref. No.: JC03014
Protect patients, nurses and doctors
from bacteria or dust
Material: 22gsm, 38gsm PP non-woven
Sizes: L-115cmx135cm,
XL-125cmx145cm,
XXL-130cmx150cm
Color: blue, yellow, pink



Non-woven Round Cap

Ref. No.: JC03015
Material: PP non-woven, single elastic
Size: 18", 21", 24"
GSM Weight: 10g
Color: blue, white
Used in medication, light industries,
food processing etc.



Surgical Gown

Ref. No.: JC03022
Protect doctors during surgical processes.
Material: PP, SMS, PP+PE, Spunlace,
Microporous, Reinforced
GSM Weight: 40g, 45g
Size: S: 130cmx115cm
M: 140cmx115cm
L: 145cmx120cm
XL: 150cmx130cm



Non-woven Strip Cap (Mob Cap)

Ref. No.: JC0301501
Material: PP non-woven, single elastic
Size: 20", 21"
GSM Weight: 10g
Color: blue, white
Used in medication, light industries,
food processing etc.



Non-woven Bed Sheet

Ref. No.: JC03011
Material: 30-40 gsm PP non-woven,
Size: 60cmx180cm, 100cmx200cm
Color: blue, white



Non-woven Doctor Cap

Ref. No.: JC03016
Material: 30 gsm PP non-woven
Size: 62cmx13cm
Forehead size: 11cm or 21cm
Color: blue, green
*with tie



Non-woven Pillow Cover

Ref. No.: JC03012
Material: 30-40 gsm PP non-woven,
Size: 60cmx180cm, 100cmx200cm
Color: blue, white



Shoe Cover

Ref. No.: JC03017
Material: 20-40 gsm PP non-woven,
PE 2.5gsm
Size: 39cmx15cm
Color: blue



BANDAGE SERIES

<http://www.joycrown.com>



Elastic Adhesive Bandage (EAB Bandage)

Ref. No.: JC04001

Material: 100% cotton elastic
• Strong Adhesion
• Easy unwind
• Center line for even overlap
Length: 4.5m
Width: 5cm, 7.5cm, 10cm



Elastic Crepe Bandage

Ref. No.: JC04005

Material: 80% cotton, 20% spandex
Weight: 65g/m², 75g/m², 85g/m²
Length (stretched): 4m, 4.5m
Width: 5cm, 7.5cm, 10cm, 15cm
Color: raw white, bleached
Packing: individual



EAB Raw Material

Ref. No.: JC04002

Raw Material for EAB Bandages,
Jumbo roll
Weight: 140g/m²
Size: 5cm x 400m
7.5cm x 400m
10cm x 400m



Elastic Crepe bandage with Red/Blue Thread

Ref. No.: JC04006

Material: 80% cotton, 20% spandex
Weight: 65g/m², 75g/m², 85g/m²
Length (stretched): 4m, 4.5m
Width: 5cm, 7.5cm, 10cm, 15cm



Cotton Stretch Bandage

Ref. No.: JC04003

Material: pure cotton
Warp threads: >59 tex
Weight: 108g/m²
Length (stretched): 4.5m
Width: 5cm, 7.5cm, 10cm, 15cm



PBT Bandage (Conforming Bandage)

Ref. No.: JC04007

Material: 55% viscose, 45% cotton
with woven edges
Weight: 26g/m², 30g/m²
Length: 4m, 4.5m, 5m
Width: 5cm, 7.5cm, 10cm, 15cm
Packing: individual cellophane



Wool Crepe Bandage

Ref. No.: JC04004

Material: 33% wool, 67% cotton
Weight: >115g/m²
Threads per 10 cm Warp:
cotton - 54 to 66
wool - 108 to 132
Length(stretched): 4.5m
Width: 5cm, 7.5cm, 10cm, 15cm



Side-Woven Bandage (Boxing Hand Wrap)

Ref. No.: JC04008

Material: 100% cotton,
40's, 32's, 21's
Mesh: 19x15, 30x20, 30x28,
40x20, 28x24
Length: 4m, 5m, 6m, 9m, 10m
Width: 4cm, 5cm, 6cm, 7.5cm,
8cm, 10cm, 15cm, 20cm



Sterile Cotton Crepe Bandage

Ref. No.: JC0400301

Material: 100% cotton, 108g/m²,
28's, S/Z twist,
Flexibility: not less than 200%
Length: 2.5cm, 5cm, 7.5cm,
10cm, 20cm, 30cm
Width: 4.5m
Packing: individual sterile pack



Fast Dressing Bandage (First Aid Bandage)

Ref. No.: JC04009

Description: gauze / PBT bandage
with cotton / non-adherent
pad sewn inside
Length: 1m, 3m, 5m
Width: 6cm, 8cm, 10cm, 15cm
Packing: individual sterile packing





BANDAGE SERIES

<http://www.joycrown.com>



High Elastic Bandage

Ref. No.: JC040101

Material: polyester and spandex
Weight: 70g/m², 90g/m²
Length (stretched): 4.5m
Width: 5cm, 7.5cm, 10cm, 15cm
Color: skin color, bleached, blue



P.O.P Slabs

Ref. No.: JC04012

Weight: 400-420g/m²
Setting time: 2-5 min., 4-8min.
Size: 10cmx80m-4ply
15cmx80m-4ply
20cmx80m-4ply



Purflid High Elastic Bandage

Ref. No.: JC040102

Material: 52% cotton, 20% polyester 28% fiber
Weight: 80g/m², 90g/m², 100g/m²
Length(stretched): 4m, 4.5m, 5m
Width: 5cm, 7.5cm, 10cm, 15cm



Cast Padding (Orthopedic Padding)

Ref. No.: JC04013

Supporting material for surgery, swadding of P.O.P bandages.
Material: synthetic /100% cotton.
Weight: 90g/m²
Length: 2.7m
Width: 2", 3", 4", 5", 6", 8"



High Elastic Bandage

Ref. No.: JC040103

Material: 75% cotton, 25% spandex
Weight: 100g/m², 120g/m²
Length(stretched): 4m, 4.5m, 5m
Width: 5cm, 7.5cm, 10cm, 15cm
*with or without clips
*elastic band clips or metal clips



Orthopedic Casting Tape

Ref. No.: JC04014

Material: fiberglass / polyester
Weight: 230g/m²
Length: 3.6m
Width: 5cm, 7.5cm, 10cm, 12.5cm, 15cm
Packing: Individual



Skin Color Elastic Bandage

Ref. No.: JC040104

Description: high elastic bandage, plain, skin color
Material: 55% bleached cotton, 45% spandex
Length: 4m, 4.5m, 5m
Width: 5cm, 7.5cm, 10cm, 15cm



Elastic Tubular Mesh Bandage

Ref. No.: JC04015

Material: nylon 45%, spandex 55%
Size: available in different sizes
• Comfortable stretching in all direction
• No skin allergy
• Anti-slipping and freedom of movement



Plaster Of Paris Bandage (P.O.P Bandage)

Ref. No.: JC04011

Weight: 400-420g/m²
Setting time: 2-5 min., 4-8min.
Length: 2.7m, 3m, 4.8m
Width: 5cm, 7.5cm, 10cm, 15cm, 20cm



Elastic Cohesive Bandage

Ref. No.: JC04017

Material: elastic fiber and cotton
Color: various color
Weight: 30g/m², 34g/m², 78g/m²
Length: 4m, 4.5m, 5m
Width: 2.5cm, 5cm, 7.5cm, 10cm, 15cm, 20cm
Packing: individual cellophane



SURGICAL TAPES & DRESSINGS

<http://www.joycrown.com>



Zinc Oxide Adhesive Plaster (Iron Tin)

Ref. No.: JC05001
Material: fabric coated with zinc oxide adhesive. Highly breathable
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm
Length: 5m, 10m, 4yds, 5yds, 10yds
Package: 12 rolls/box
*with or without spool and cover
*zigzag edge available



Waterproof Surgical Tape

Ref. No.: JC05008
Material: waterproof fabric, hot-melt adhesive
Latex free and hypoallergenic
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm
Length: 5m, 10m, 5yds, 10yds
Package: 12 rolls/box
*with or without spool and cover
*zigzag edge available



Zinc Oxide Adhesive Plaster (Plastic Can)

Ref. No.: JC05002
Material: fabric coated with zinc oxide adhesive. Highly Breathable.
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm
Length: 5m, 10m, 4yds, 5yds, 10yds
Package: 12 rolls/box
*with or without spool and cover
*zigzag edge available



Fixed Dressing Plaster

Ref. No.: JC05007
Bulk packing, can cut into any size when use.
Material: non-woven, hot-melt adhesive
Size: 2.5cmx10m, 5cmx10m, 10cmx10m, 15cmx10m, 20cmx10m, 30cmx10m



Transparent PE Tape

Ref. No.: JC05003
Material: microporous perforated PE
Latex free and hypoallergenic
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm
Length: 5m, 10m, 4yds, 5yds, 10yds
Package: 12 rolls/box
*with or without spool and cover
*zigzag edge available



Band Aid

Ref. No.: 05008
Material: cotton fabric with highly absorbent cushion pad
Size: 72x19mm, 70x18mm, 60x18mm,
Diameter: 22mm (round shape)
Packing: various packing available
* Custom shape, size, material, color, printing available upon request



Air-permeable Surgical Tape (Paper Tape)

Ref. No.: JC05004
Material: microporous non-woven fabric, hot-melt adhesive coated
Latex free and hypoallergenic
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm
Length: 5m, 10m, 4yds, 5yds, 10yds
Package: 12 rolls/box
*zigzag edge available



Waterproof Band Aid

Ref. No.: JC05009
Material: PE with highly absorbent cushion pad
Size: 56x19 mm, 76x19 mm, 72 x 19 mm,
other sizes are available
Packing: 100pcs/box



Silk Surgical Tape

Ref. No.: JC05005
Material: acetate silk, adhesive
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm
Length: 5m, 10m, 4yds, 5yds, 10yds
Package: 12 rolls/box
*zigzag edge available



Band Aid

Ref. No.: JC05010
Material: cotton fabric or PVC with highly absorbent cushion pad
Size: 60x100mm, 60x400mm, 60x500mm
Packing: 100pcs/box



SURGICAL TAPES & DRESSINGS

<http://www.joycrown.com>



I.V. (Intro Venous) Dressing
 Ref. No.: JC05011
 Fixing and caring for I.V. cannular and all kinds of infusion lines.
 Material: microporous non-woven fabric, hypoallergenic adhesive, absorbent pad,
 Size: 6cmx7cm, 6cmx8cm, 7cmx8.5cm, 8cmx10cm
 Packing: individual sterile pack



PU Incise Drape
 Ref. No.: JC05016
 Use for bone abdomen operation, chest and heart surgical operation, plastic surgery, and other operations
 Material: PU/PE
 Size: 15cmx28cm, 15cmx20cm, 45cmx55cm, 90cmx90cm
 Packing: individual sterile pack



Kinesiology Tape
 Ref. No.: JC05018
 - Supporting bandage for sore muscles;
 - Assist in lymphatic drainage;
 - Activate endogenous analgesic system;
 - Correct joint problem
 Size: 2.5cmx5m, 5cmx5cm



Sports Medicine Adhesive Tape
 Ref. No.: JC05017
 Material: fabric, zinc oxide adhesive or hot-melt adhesive coated
 Color: bleached white, black
 Size: 1"x10 yards (2.5cm x 9.1m)
 2"x15 yards (5.0cm x 13.7m)
 *zigzag edge available



Self Adhesive Wound Dressing
 Ref. No.: JC05013
 Cover and protect wounds after operation.
 Material: non-woven fabric, anti-allergenic adhesive, absorbent pad
 Size: 7x5cm, 8x9cm, 9x10cm, 9x15cm, 9x20cm, 9x25cm, 9x30cm
 Packing: individual sterile pack



Calcium Alginate Wound Dressing
 Ref. No.: JC05023
 Made of alginate which is brought from natural seaweeds.
 Use on skin, trauma and other refractory wound.
 Easy to use, good air permeability, excellent biocompatibility.
 Can be absorbed by human body.



Transparent Dressing
 Ref. No.: JC05014
 Waterproof, impermeable to bacteria, hypoallergenic
 Material: PE film, hypoallergenic adhesive
 Size: 4cmx5cm, 6cmx7cm, 7cmx11cm, 10cmx12cm, 15cmx20cm



Wound Closure Strip
 Ref. No.: JC05022
 Material: non woven fabric, coated with hypoallergenic adhesive.
 Provides wound support and assists in increasing the tensile strength of wound.
 Size: 38x6mm, 75x3mm, 75x6mm, 100x6mm, 100x12mm
 Packing: 3/5/6 strips/pack



Hydrocolloid Dressing
 Ref. No.: JC05015
 Used by special lipidic hydrocolloid technology, spread hydrocolloid on high molecular polymer and acetic acid non-woven fabric covered by cosmoil oil.
 Use for: skin trauma, burn, scald
 Size: 5cmx5cm, 7.5cmx7.5cm



Perforated Z.O.P
 Ref. No.: JC05021
 Material: skin color fabric coated with zinc oxide adhesive, perforated.
 Highly breathable.
 Color: skin color
 Size: 5cmx5m, 10cmx5m, 18cmx5m
 Packing: 1 roll/plastic bag or box.



SYRINGES

<http://www.joycrown.com>



Disposable Syringes with needle (Long blister packing)

Ref. No.: JC06001
Material: medical polymer, non-toxic, pyrogen free
Size: 1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 50ml
Needle: 18-26G, 16G, 27G,
Packing: individual blister pack
*sterile



Hypodermic Needles

Ref. No.: JC06003
Latex free. Sterile, Non-toxic, pyrogenic free
18/8 stainless steel
Size: 18-26G, 16G, 27G, 29G
Packing: individual blister pack
*sterile



Disposable Syringes with needle (Short blister packing)

Ref. No.: JC0600101
Material: medical polymer, non-toxic, pyrogen free
Size: 1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 50ml
Needle: 18-26G, 16G, 27G,
Packing: individual blister pack
*sterile



Disposable Infusion Set

Ref. No.: JC06004
Material: medical polymer, non-toxic, pyrogen free
Parts: Non-vented plastic spike, Inject-molded chamber, PVC tubing, PP regulator, PP luer slip connector.
Packing: individual blister pack or poly pack
*sterile



Disposable Insulin Syringes

Ref. No.: JC06002
Material: medical polymer, non-toxic, pyrogen free
Main Sizes: 0.5ml, 1ml, 3ml
Needle: 18-26G, 16G, 27G, 29G
Color: Orange
Packing: individual blister pack
*sterile



I.V. Catheter

Ref. No.: JC06005
Material: High quality needle, PTFE & FEP & PUR tube
Size: 14g, 16g, 18g, 20g, 22g, 24g
* with or without wings
* with or without injection valve
Packing: individual blister pack
*sterile

URINE BAGS



Urine Bags

Ref. No.: JC09001
Material: Non toxic, medical grade PVC
Size: 1000ml, 2000ml
Packing: individual sterile pack.
Outlet: T outlet or pull-push outlet



Pediatric Urine Collector

Ref. No.: JC09004
Material: Non toxic, medical grade PVC
Size: 100ML
Packing: individual sterile pack.



TEST KITS

<http://www.joycrown.com>



HCG Pregnancy Test Strips

Ref. No.: JC07001

Designed for professional and home users.
99% accurate results in under 5 minutes.
Packing: individual pack with a moisture proof pouch.



LH Ovulation Test Strip

Ref. No.: JC07004

The LH Ovulation Test is specifically designed to detect your LH surge - the time when you are likely to ovulate. If you receive a positive on an LH test, you are in your most fertile phase of your menstrual cycle.
Packing: individual pack with a moisture proof pouch



HCG Pregnancy Test Cassette

Ref. No.: JC07002

HCG pregnancy test in cassette format is highly accurate. This format cassette pregnancy test is used by hospitals and clinics.
Packing: Individual pack with a moisture proof pouch.



LH Ovulation Test Cassette

Ref. No.: JC07005

Ovulation urine test cassette is mostly used by professionals or home users because it is better protected from mistakes compare to dip strips. Simply place several urine drops into cassette and wait 5 minutes for results.
Packing: individual pack with a moisture proof pouch.



Pregnancy Test Midstream

Ref. No.: JC07003

This format is primarily design for home use. These tests are placed directly in the urine stream, there is no need to collect urine in a cup. The midstream pregnancy test is over 99.9% accurate.
Packing: individual pack with a moisture proof pouch.



LH Ovulation Test Midstream

Ref. No.: JC07006

The LH Ovulation Test Midstream will help you find the best time when you are most able to conceive (become pregnant). It can help you plan your pregnancy.
Packing: individual pack with a moisture proof pouch.

THERMOMETERS



Clinical Thermometer

Ref. No.: JC08001

To test body temperature
Temperature range:
35 - 42 degrees Centigrade
94-108 Fahrenheit degrees
Accuracy: Celsius 0.1
Fahrenheit 0.2
Packing: individual plastic box



Digital Thermometer

Ref. No.: JC08002

Thermometers offer fast, safe and reliable temperature readings. An audible beep signals the completed measuring process once peak temperatures are reached.
Head type: soft or hard
Packing: individual plastic box



FIRST AID KITS

<http://www.joycrown.com>



Surgical Dressing Kits (Pack)

Ref. No. JC10001
Each pack contains:
1 x Tray
1 x Fine Forceps 14cm
1 x Plastic Tweezers
5 x Compress 20 x 20cm folded
5 x Lint (Gauze Balls or cotton)
1 x Surgical Drape 50 x 60cm



Surgical Dressing Kits (Pack)

Ref. No. JC10002
Each pack contains:
1 x Tray
2 x Fine Forceps 14cm
5 x Compress 20 x 20cm folded
5 x Lint (Gauze Balls or cotton)



Sterile Dressing Pack

Ref. No.: JC10003
Specification No. 10 & 35
Each Pack Contains:
1 x Absorbent Tissue Pad
4 x Absorbent Swabs, 10cmx10cm
4 x Large Absorbent Cotton Balls
1 x Absorbent Paper Towel
1 x Water Repellent Paper Sterile Field



First Aid Kit

Ref. No.: JC10004
Each pack contains:
1 x Triangle bandage
1 x Zinc Oxide Adhesive Plaster, 1.25x5m
2 x Adhesive plaster, 6x10cm
4 x Sterile gauze folded, 7x10cm
1 x Sterile cutting gauze, 6x6cm/100pcs
4 x Sterile fasten bandage, 6cmx5m
4 x Sterile fasten bandage, 10cmx5m
1 x Scissors



First Aid Kit

Ref. No.: JC10005
Each pack contains:
1 x Adhesive tape
5 x Bum Dressing
1 x Tweezers
1 x Emergency Blanket
1 x Scissors
1 x Crepe Bandage
2 x Triangular Bandage
4 x Disposable Gloves
8 x Adhesive Plaster, 10cm x 5cm
1 x ABS Plastic Box with OEM label or printing



First Aid Kit

Ref. No. JC10006
Each pack contains:
10 x Adhesive Plaster
10 x Cotton Swab
1 x Scissors
1 x Non Woven Tape
2 x Gauze Bandage: 5CM
2 x Gauze Bandage: 7.5CM
2 x Soap Wipes: 12X22CM



First Aid Kit

Ref. No. JC10007
Each pack contains:
5 x Adhesive Bandage
4 x Non Woven Pad: 5cm x 5cm/4ply
2 x Alcohol Pad
2 x BZK Antiseptic Swab
2 x PBT Bandages
1 x PVC Glove: L Size
1 x Adhesive Tape: 1.25cm x 1m
1 x Scissors
10 x Safety Pin: 1 x Tweezers
1 x Emergency Blanket
2 x First Aid Bandage: 1 x Eye Pad



Skin Traction Kit

Ref. No.: JC10008
The extension strapping in this kit is coated with a skin-friendly adhesive, to which soft foam has been bonded to protect the foot.
Size: adult, child
Packing: individual printing box



Bag for First Aid Kits

Ref. No.: JC10009
Material: 210D, 420D, 600D nylon
Size: any size and design available as request



Emergency Blanket (Rescue Sheet)

Ref. No.: JC10010
Keep body heat, waterproof, for maximum safety.
Material: PET, 12micron thick
Size: 140cmx210cm; 160cmx210cm
Color: silver or gold
Packing: individual PP pack



ACCESSORIES

<http://www.joycrown.com>



Mouth to Mouth Quick Saver (Rescue Mask)
Ref. No.: JC10011

Material: plastic, with one-way valve
Size: 20cmx20cm
Packing: individual zip bag pack



Surgical Scissors
Ref. No.: OT006

Bandage Scissors
Iris Curved Scissors
Size: total length 18cm or 4"



Hot & Ice Bag
Ref. No.: JC10012

For warm and effective relief from muscular and joint aches and pains. For an effective cooling treatment, medical and sport
Size: 12 x 29cm
10cm x 13cm



Plastic & Metal Forceps
Ref. No.: JCOT007

Material: plastic, metal
Size: 10.8cm, 11.5cm, 12.5cm, 14cm



Cervical Collar
Ref. No.: JC11001

Material: high density closed cell foam (100D, 40D, 60D), with 100% cotton woven covering
Color: white, black
Size:
Small (S): length 44cm; width 9cm
Medium (M): length 48cm; width 10cm
Large (L): length 53cm; width 11cm
Packing: individual pack with PP bag and box



Plastic Medicine Cups
Ref. No.: JC13002

Material: PP
Weight: 1g, 1.2g, 1.4g or others
Size: 30ml, 60ml, 120ml



Cervical Cuff
Ref. No.: JC11002

Material: foam, pure cotton
Size: 5cm x 6m; 7.5cm x 6m
Package: 2 rolls into one box within 12pcs fastener



PVC Examination Gloves
Ref. No.: JCOT009

Material: PVC
Latex free, non allergic
Thickness: min.0.05
Tensile Stretch: >300
Size: S(85cm), M(95cm), L(105cm), XL(115cm)
*powdered or powder free



Neck Band
Ref. No.: JC11003

Tube holder with adherent fastener
Size: 30cmx3cm
* other size available as request



Surgical Gloves
Ref. No.: JCOT010

Material: medical grade latex
Size: 6', 6.5', 7', 7.5', 8', 8.5', 9'
*powdered or powder free
*sterile

SHANGHAI JOY CROWN INDUSTRY CO., LTD.

**Add.: Rm. 220A, Bld.11-C, No.613 E'Shan Road, Pudong New Area
Shanghai 200127, P.R. China**

Tel.: +86-21-58751019

Fax.: +86-21-50907607

Email: sales@joycrown.com

[Http://www.joycrown.com](http://www.joycrown.com)